

# Experiment 7


## Flame Photometry

# INTRODUCTION

Atomic spectroscopy is based on the absorption, emission, or fluorescence process of light by atoms or elementary ions. Information for atomic scale is obtained in two regions of the electromagnetic radiation (EMR) spectrum. These regions are UV/VIS and the X-ray.


**Atomic spectroscopy can be divided as:**

1. Absorption spectroscopy (UV-Vis)
2. Emission spectroscopy (Flame photometry)
3. Luminescence / fluorescence spectroscopy




When sample interacts with light, absorption process occurs. Ground state electrons of the sample atom tend to move to the excited states with the energy of absorbed light. This process can also be called **excitation**.

- Apart from light, **heat** can cause excitation. Since excited state is unstable, electrons want to return back to the ground state. When an excited electron turns back to its ground state a radiation is emitted that is equal to the energy difference between excited and ground states. The emitted light is monochromatic, and it has the same wavelength as the light absorbed in the excitation process.



- Depending on the excitation technique, absorbed or emitted light is measured. If excitation source is flame, emitted radiation is measured. On the other hand, absorption is measured when lamp is used for excitation. Both are directly proportional with the number of atoms in the sample.



# DEFINITION

**Flame photometry** or **flame emission spectroscopy** is an atomic emission technique used for the determination of elements which can be easily excited to higher energy levels at flame temperature such as **alkali (Cations) and alkaline earth metals**. These metals are Na, K, Ca, Ba, Li.



LITHIUM



SODIUM



POTASSIUM



CALCIUM



BARIUM

This method is based upon the measurement of **intensity of radiation emitted**, in the **visible region**, when a metal atom is introduced into a flame. The wavelength of the radiation (or the color), emitted tells us what the element is (**qualitative**), and the intensity of the radiation tells us how much of the element is present (**quantitative**).

# Flame

*A flame can be described as a steady state gas phase reaction which takes place with emission of light. It is produced by burning a mixture of **fuel** and **oxidant** in a burner. In flame photometry a variety of **fuels** can be used. Oxidant can generally be air, oxygen or nitrous oxide ( $\text{N}_2\text{O}$ ). The flame temperature depends on fuel to oxidant ratio.*



*❖ The most commonly employed flames in flame photometry can be grouped into two types:*

1. Flames in which the fuel and oxidant as air or oxygen are well mixed before combustion, these are called **pre-mix** or **laminar flames** as they exhibit laminar flow.
2. Flames in which the fuel gas and the oxidant are first mixed in the flame itself. They are called **unpremix** or **turbulent flames** since they exhibit turbulence

Flame is formed by two components: **fuel** and **oxidant**.

**Temperature of the flame changes** depending on the fuel and oxidant types and their proportions.

In flame photometer generally **natural gas is used as a fuel** and **air is the oxidant**. Table 1 lists the different types of fuel, oxidant, and the temperature of the flame.


| FUEL        | OXIDANT       | TEMPERATURE, °C   |
|-------------|---------------|-------------------|
| Natural Gas | Air           | <b>1700</b> -1900 |
| Natural Gas | Oxygen        | 2700-2800         |
| Hydrogen    | Air           | 2000-2100         |
| Hydrogen    | Oxygen        | 2550-2700         |
| Acetylene   | Air           | 2100-2400         |
| Acetylene   | Oxygen        | 3050- <b>3150</b> |
| Acetylene   | Nitrous Oxide | 2600-2800         |



# Flame Structure

Flame consists of four important regions that their appearance and relative sizes changes with the fuel-oxidant ratio. These are:



- i. **preheating zone** where the combustion mixture is heated to the ignition temperature by thermal conduction from the primary reaction zone.
- ii. The **primary combustion zone** of the flame is **blue** in color. In this region, the concentration of ions and free radicals is very high and there is no thermal equilibrium. Therefore, it is not used in flame spectroscopy.
- iii. The **interconal region** Where the maximum temperature is achieved just above the tip of the inner zone. It is rich in free atoms and is the most widely used region for flame spectroscopy.
- iv. The **outer cone** (secondary reaction zone) - In this zone, the products of the combustion processes are burnt to stable molecular species by the surrounding air.






# TECHNIQUE

In this technique, first aerosols are formed from sample solution by a jet of compressed gas. This process is called **nebulization**. Then the flow of the gas carries the aerosols into a flame where atomization takes place. **Atomization** is the conversion of sample aerosols into an atomic vapor by flame.







When a sample is atomized, a substantial fraction of the metallic constituents is reduced to gaseous atoms. Depending on the temperature of the flame a certain fraction of these atoms is ionized. Then, electrons of the formed atoms are excited to upper state. Light is emitted at characteristic wavelengths for each metal as the electron returns to the ground state.

# Flame Photometer Systems


In flame photometer there are three fundamental systems which are **emission**,  **$\lambda$ -selection**, and **recording systems**.



1) **Emission System:** This consists of the flame, which is the source of emission.



| Element       | Emitted wavelength | Flame colour                                                                                     |
|---------------|--------------------|--------------------------------------------------------------------------------------------------|
| Potassium (K) | 766 nm             | Violet        |
| Lithium (Li)  | 670 nm             | Red           |
| Calcium (Ca)  | 622 nm             | Orange        |
| Sodium (Na)   | 589 nm             | Yellow       |
| Barium (Ba)   | 554 nm             | Lime green  |



**2)  $\lambda$ -Selection System:** This includes the whole optical system of wavelength selection. In flame photometer the wavelength selector is a **filter**.


- The radiation emitted by the excited atoms is selected by using a filter which **transmits an emission line** of one of the elements **while absorbing the others**.
- **There are two types of filters.** These are **absorption** and **interference** filters. Absorption filters are restricted to visible region of the spectrum, but interference filters are used in UV, VIS and IR regions of the spectrum.

- **Absorption filters** are less expensive than the interference filters and they have been widely used for band selection in the visible region. These filters function by absorbing certain portions of the spectrum and transmitting the band of wavelengths belonging to the analyte element. The most common type consists of colored glasses.
- **Interference filters** rely on optical interference to provide relatively narrow bands of radiation. They consist of a transparent dielectric layer ( $\text{CaF}_2$  or  $\text{MgF}_2$ ) that occupies the space between two semi-transparent metallic films. This array is sandwiched between two plates of glass.



3) *Recording System*: This part consists of all the means of detection (*phototubes* or *photomultiplier tubes*), the *electronic devices of amplifying* and *electrical apparatus* for *measuring and direct recording*.


## Flame Photometry process



# The BWB XP Flame Photometer Instrument


The instrument that is used in this experiment is a **BWB XP model**, a five-channel digital flame photometer which is a low temperature, designed for the routine determinations of Na, K, Ca, Li, Ba. It is a direct reading digital instrument designed for use in clinical, industrial, and educational applications





# Calibration Curve

In flame photometry emitted light intensity from the flame is directly proportional to the concentration of the species being aspirated. The graph below shows that the direct relationship between **the emission** and **concentration** is true only at relatively low concentrations of mg/L level (up to 50 mg/L) → (ppm).



- ❖ If the samples being analyzed lie on the linear part of the curve, then user can take direct concentration readings from the digital display. However, if the concentration of sample is above the linear part, then user must dilute sample so that it lies on the linear part of the curve.
- ❖ A calibration curve is obtained by using standard solutions containing known concentrations of the elements to be determined.
- ❖ The concentration range covered by the calibration curve depends on the expected concentration so that the sample readings fall somewhere inside the calibration curve.
- ❖ Once the calibration curve has been plotted, the scale reading for the sample solution is compared with the curve to find the concentration.
- ❖ It is important to emphasize that each element has its own characteristic curve and separate curves must be constructed.



A laboratory setting featuring a blue flame photometer on the left with a digital display showing '20.0ppm', '18.0ppm', '3.0ppm', and '12.0ppm'. To its right are several glass vessels: a small Erlenmeyer flask, a round-bottom flask, a large Erlenmeyer flask, and a beaker, all containing dark liquids. In the background, a rack of test tubes and a computer monitor are visible. A semi-transparent blue banner with the text 'Practical Part' is centered over the image.

# Practical Part

# Objective


- ❖ Determination of **Na** and **K** assay in oral rehydration solution (ORS) as a pharmaceutical product using flame-photometer.
- ❖ Get familiar with the BWB Flame -photometer instrument and its usage.



Oral rehydration therapy is a type of fluid replacement used to prevent and treat dehydration, especially due to diarrhea. It involves drinking water with modest amounts of sugar and salts, specifically sodium and potassium.

➤ Unlike other fluids, the ratio of the ingredients in an ORS matches what the body needs to recover from a diarrheal illness.

|                                                                                                                                      |         |                      |
|--------------------------------------------------------------------------------------------------------------------------------------|---------|----------------------|
| <b>Composition:</b>                                                                                                                  |         |                      |
| Each sachet contains                                                                                                                 |         |                      |
| Sodium Chloride                                                                                                                      | 0.520 g | Na+                  |
| Potassium Chloride                                                                                                                   | 0.300 g | K+                   |
| Tri Sodium Citrate Dihydrate                                                                                                         | 0.580 g | Cl-                  |
| Anhydrous Dextrose                                                                                                                   | 2.700 g | Citrate              |
|                                                                                                                                      |         | Anhydrous Dextrose   |
|                                                                                                                                      |         | Total Osmolarity     |
|                                                                                                                                      |         | 245                  |
|                                                                                                                                      |         | m.mole concentration |
|                                                                                                                                      |         | 75 mmol / liter      |
|                                                                                                                                      |         | 20 mmol / liter      |
|                                                                                                                                      |         | 65 mmol / liter      |
|                                                                                                                                      |         | 10 mmol / liter      |
|                                                                                                                                      |         | 75 mmol / liter      |
| <b>Dosage :</b> As directed by the physician.                                                                                        |         |                      |
| Store below 25°C, protect from light and moisture.                                                                                   |         |                      |
| <b>Direction for use :</b>                                                                                                           |         |                      |
| * Cut open the sachet , empty the powder in 200 ml of drinking water and mix well. Do not boil the solution.                         |         |                      |
| * Orally feed the solution to the patient in small amounts as per the need. In case of vomiting, wait for 10 minutes and give again. |         |                      |
| * Refer the patient to a doctor in case of persistent vomiting and or continuous diarrhea.                                           |         |                      |
| * In between the feeds, supplement by other liquids (breast milk in case of infants), fruit juice and light food.                    |         |                      |
| * The reconstituted solution should be used immediately within 24 hours and should not be stored.                                    |         |                      |
| NCCJ 06                                                                                                                              |         |                      |





# Glassware

Each group will be supplied by the following glassware:



# Procedure

**Step 1:** Preparing the ORS solution by dissolving the powder in one sachet in 250ml water. Then a dilution of 5ml to 100ml will be carried out.




## Step 2: Preparing Standard Calibration curves


Preparing calibration curve for (Na and K) using NaCl and KCl salts. At first two stock solutions (250ml) with the concentration **200ppm** from both metals. Then carrying serial dilutions to obtain 4 standards with the concentrations (**100ppm**, **80ppm**, **40ppm**, and **20ppm**).



Flame photometer Na Standard graph



Flamephotometer potassium K standard graph




Weigh a certain amount of salt



200ppm


Dilute 25 ml from the stock solution to 50 ml D.W



100ppm

S1


Dilute 10 ml from the stock solution to 25 ml D.W



80ppm

S2

Dilute 10 ml from the stock solution to 50 ml D.W




40ppm

S3

Up to 250 ml D.W



Stock solution




Step 3: Take the standards to the flame-photometer to be measured



Serial Dilution  
(2 steps dilution)

Then Dilute 25 ml from standard 3 (S3) to 50 ml D.W



20ppm

S4

# How to prepare a stock solution in ppm?

**Example:** Prepare 250ml solution of Ca with a concentration of 100ppm

$$\text{ppm} = \text{mg/L}$$

$$100\text{ppm} = 100\text{mg/L} \rightarrow (250\text{ml} = 0.25\text{L})$$

$$100\text{mg} \rightarrow 1\text{ L}$$

$$??\text{Ca} \rightarrow 0.25\text{L}$$

$$= 25\text{mg} (0.025\text{g}) \text{ of Ca per } 250\text{ml}$$

How much should we weight of  $\text{CaCl}_2$  ??



## **\*\*Equivalent weight\*\***

0.025g Ca in 250ml → 100ppm Ca



1 mol of Ca = 1 mol of CaCl<sub>2</sub>

$$\text{Mol} = \text{weight} / \text{MWT}$$

$$0.025\text{g} / 40.078 = ??\text{g} / 110.98$$

$$0.00062 = ??\text{g} / 110.98$$

$$\text{Weight of CaCl}_2 = 0.06881 \text{ g} = 68.81 \text{ mg}$$

❖ Molecular weight for  
(Ca= 40.078 g/mol)  
(Cl= 35.453 g/mol)  
(CaCl<sub>2</sub>= 110.98 g/mol)



# Sodium from NaCl

250ml of  $\text{Na}^+$  in a  
concentration of  
(200 ppm)

ppm = mg/L

$\text{Na}^+ ??$


Equivalent Wt.

$\text{NaCl} ??$



❖ Molecular weight for  
(Na= 22.990 g/mol)  
(Cl= 35.453 g/mol)  
(NaCl= 58.44 g/mol)

200 ppm of  $\text{Na}^+$



# Potassium from KCl

250ml of  $K^+$  in a  
concentration of  
(200 ppm)

ppm = mg/L

$K^+ ??$


Equivalent Wt.

$KCl ??$





❖ Molecular weight for  
(K= 39.0983 g/mol)  
(Cl= 35.453 g/mol)  
(KCl= 74.55 g/mol)

200 ppm of  $K^+$




# How to Obtain the accurate volume (liquid meniscus)

**Remember**



- ✓ After finishing measurements, draw calibration curve then calculate the concentration for your Na and K metals in the sample(D.f).
- ✓ Determine the ORS sample Assay for both (Na and K) and discuss wither you are with or against its claim.
  - To determine the % assay, use the following formula:

$$\% \text{ Assay} = \frac{\text{Actual Concentration}}{\text{Theoretical Concentration}} \times 100\%$$



| m.mole concentration |                 |
|----------------------|-----------------|
| Na+                  | 75 mmol / liter |
| K+                   | 20 mmol / liter |
| Cl-                  | 65 mmol / liter |
| Citrate              | 10 mmol / liter |
| Anhydrous Dextrose   | 75 mmol / liter |
| Total Osmolarity     | 245             |



GOOD  
LUCK!

*MSc. Farah Hudaib*