


Quantitative and Qualitative Analysis of Ascorbic Acid

water soluble vitamine antiaging ك تبه له مورد دون الله ما مد Ascorbic acid (well-known as vitamin C) is products a naturally occurring weak acid, with antioxidant properties and an important component of a healthy diet. is It is needed by the body to help wounds heal, to Vite 205 Fix 200 Supplements 1 * enhance the absorption of iron from plant Comula cina formula cina foods, and to support the immune system. - important protein in cartilage, tendens, shin, blood vessels, It is essential in collagen synthesis and works as an antioxidant to protect cause damage to skin your cells against free radicals, which may play a role in heart disease, which may play a role in heart disease, antiaging cancer, and other diseases. The history of Vitamin C revolves around the antiwrinkles who seed Vit Crab Sedi history of the human disease scurvy. Its symptoms include exhaustion, agent massive hemorrhaging of flesh and gums, general weakness, and diarrhea.

cais, & antioxidant is vitamine (1 aily apoil) to protect your cells against free radicals grand of con the skin Esteris neutrolization un de ascorbic acid sos só aging much الانكرونات الله وحاجها ملا عاج الجد Cancer of heart diseasement is so مرجدا وعجدا انه VitC ما اله فوات ما ناحية وقائمة رما جمة للعراق معلومة المحمور من المحمور المجامد على المحمور المعامد) المحمور المعامد) المحمور المعامد) xfree radicals voisse وإلى Vitamiae Cسام بعطيها هاد الالحكرن إلى chemicale ho, pullation he entre entre la source de la constante UV radiation les means, الحجامية . وفي فيتامينات موجودة بالصيلات ماعة الشم (col gos (anti oxidants) (parly) مهور کیثر فی واحکات الش الفائمة المطبحة الكاسة grand VitC as is spir Fe suppliment old is * non absorbable form. Absorbably of cits sleek as baid - Recious splend se and ferrice -> ferrous Supporting innate d'all Mix immunity system & boost & by - stimulation of reutrophile & site of in Rection improve . P phagocytosis and oxidant generation and in microbial age Killing. من الماد تسور الما من المادي والما المرب المادي ال الناصة الجعيلية مناسبة (المكتورة شرحتم).

e radicals, and other

ds heal, to ods, and to






to bran antiaging , cancer, and other diseases. The history of Vitamin C revolves around the Scarve was - with VitCrab seed, antiwrinkles history of the human disease scurvy. Its symptoms include exhaustion, agent massive hemorrhaging of flesh and gums, general weakness, and diarrhea. ا ي مهم العداد . stille siz zer development of red & blue spots on the skin Resultant death was common. Jesus Tones cos usually in your shins have skin that bruises Vitamin C is a six-carbon chain, closely related chemically to glucose. المله محودة من سلايد In 1934, Rechstein worked out a simple, inexpensive, four-step process رك شرصا كا for synthesizing ascorbic acid from glucose. This method has been الدكسورة فحدثابة Symptoms 11 used for commercial synthesis of Vitamin C. - List adoute al philosophicats of the sales the first glucose anie e Ascorbic Acid (Structure and Formula) proben donor - antioxidant by neutro lizing the free radicals in bodu * اي جرعا = زيارة عه ط hydrayl groups OR JOXIC JUSCINIE that give it a water soluble criteria OH cim Eso Ag any excess of VitC will not cause any problem become Polarity when II its easily excrehed

The history of Vitamin C revolves around the history of the human disease scurvy. Its symptoms include exhaustion, massive hemorrhaging and swollen of flesh and gums, loss of teeth, general weakness. diarrhea, development of red or blue spots on the skin, usually on your shins have skin

that bruises easily. Resultant death was common if left untreated .



الي تحته خط مو موجود بملف السلايدات عاليتمز بس هاي السلايده شرحت عليها الدكتورة و ذكرت النقطه هاي



Vitamin C can be found in various products ranging from supplements to skincare products.





* VitC has several dosage Rom.

is a sice 25i al

D) gumes

2) lotions medisjointh

3) Serum

4) additive to sunscreen J for fine lines/winkles

5) Cream

6) drops

* Vit C should be obtained from outer source

المال المال

pale skin (2)

powerful reducing agent
free radiculus = 1,751 was
Solution pails Properties of Ascorbic acid:

Naturally occurring organic compound (weak acid). - mildy a cidic solon. Obs

White to yellow crystals or powder.

Water-soluble vitamin (mildly acidic solution).

(unstable)

In the solution → rapidly oxidized (destroyed) upon exposure to heat, light, metals, and oxygen. (Cu) Cuppor mainly

In the dry state → reasonably stable to air (Not to UV light). لازا مفطة بعكان مهم وعكى سترول لا عامدهم وعكى سترول

Vitamin C occurs naturally primarily in fresh fruits and vegetables.

(ما	السفا السفا
Vitamin-C (mg/100g)	Citric acidliagis i le VII Card Gines
100 - 350	Chili peppers, sweet peppers, parsieg, and turnip greens
25 - 100	Citrus juices (oranges, lemons, etc.), tomato juice, mustard greens, spinach, brussels sprouts
10 - 25	Green beans and peas, sweet corn, asparagus, pineapple, cranberries, cucumbers, lettuce
<10	Eggs, milk, carrots, beets, cooked meat



ascorbate of regaline charge becouse of resonance system

البوتون الموتون الموت

In ascorbic acid structure, you can see the conjugated system, composed of the electrons in the double bond, hydroxyl group lone pair, and carbonyl

group. Because of resonance stabilization of the conjugate base (stable) (ascorbate), the hydroxyl group in ascorbic acid is much more acidic than typical hydroxyl groups.

المَّارِينَ عَلَيْ مِلَا لَكُورُكُ اللَّهِ الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِ اللَّهِ اللَّهِ الْمُعَالِمِي الْمُعِلَّمِ الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعَالِمِي الْمُعِلَّمِي الْمُعَالِمِي الْمُعِلَّمِي الْمُعِلَّامِي الْمُعِلَّمِي الْمُعِلَّمِي الْمُعَالِمِي الْمُعِلَّمِي الْمُعِلِمِي الْمُعَا

assorbate (stable)

titration: always occur in pairs of electrons. gain of elections decrease in oxidation number Reduction Oxidation احتر ال تأكس کسالکتر بنا = فقد الكروتات شحنة تزاد شحنة تقل ملا واصيفتع لا ١٦١١ كاللي سكيب loss of electrons Increase in oxidation number

Quantitative Determination of Ascorbic Acid: where in free radialed - l'ensil slaste free radicule sip

> A suitable method for the determination of vitamin C (ascorbic acid, C6H8O6) is a titration with potassium iodate (KIO3). In this reaction (redox titration), potassium iodate is used as a titrant and when added to an ascorbic acid solution that contains strong acid and potassium iodide (KI), the potassium iodate reacts with potassium iodide, liberating molecular iodine (I2) as shown in the reaction below:

detau ul Courter pre la [1] $\underline{\text{KIO}_3} + \underline{5\text{KI}} + 6\text{H}^+ \rightarrow 3\text{I}_2 + 6\text{K}^+ + 3\text{H}_2\text{O}$ axorbic acid most of the reduce oxidizing agent Iz desir and last on

During the titration, if the solution contains ascorbic acid, the I2 produced in equation 1 is used up in a rapid reaction with ascorbic acid (equation 2), during

@ redox RXV, delegel 4 which dehydroascorbic acid (C₆H₆O₆) and iodide ion (I-) are formed: من تنكل (1) Litration with $C_6H_8O_6 + I_2 \rightarrow C_6H_6O_6 + 2I^- + 2H^+$, was soon Jascorbic HO Ascorbic acid 2Ioxidizing agent reduction Once all the ascorbic acid has been consumed, any excess iodine (I_2) will remain in solution. This excess iodine reacts with starch, to form an intensely blue colored

emplex, indicating that the endpoint is reached. care outline of endpoint of

علام المون الرق مع المون البرعة على المون أرزى م ح اللون البرعة الون البرتقال بتجربتنا بيطلع مايل الموك الأخضر

Redox

titration

Practical Work

Materials & Chemicals Needed:

- Citric Fruit (Orange or Limon).
- Starch indicator solution: (1%).
- H2SO4 (1M).
- Iodine Solution (0.05M).
- Vitamin C tablet.
- Bendict's reagent
- NaHCO3
- AgNO3

Glassware Needed:

- Burette and stand
- 100 mL volumetric flask
- 20 mL pipette
- 250 mL Erlenmeyer flask
- Dropper.
- Graduated cylinders (10ml) & (100).
- Test tubes
- Spatula

Notes:


· Be cautious while dealing with acids they can cause burns and irritation, so avoid direct contact and keep it in the fumehood.

Part 1: Determination of Ascorbic Acid Amount (mg%) in Fruit Juice:



Procedure:

- 1. Freshly juice a citric fruit (an orange) to obtain its juice.
- 2. Take 10 mL of the fruit Juice, place in a 250 mL Erlenmeyer Flask.
- 3. Add 80 mL of distilled water and stir.
- 4. Add 10 mL of 1M H₂SO₄ (Fume hood).
- 5. Add 5 mL of starch solution (as indicator).
- 6. Start titration with 0.05M lodine solution until you reach the end point.
- 7. Record the volume of iodine solution (titrant) used to reach the end point.
- 8. Calculate the (%mg) of ascorbic acid present in your fruit sample.



> The sequence of colors is clarified in the following figure:

 Left flask: before the endpoint, the color of the solution reflects the bright orange color of fresh orange juice.

Centre flask: Once all the ascorbic acid has been oxidized, a slight excess of added iodate forms a starch-iodine complex, a green color in this case. This is the endpoint of the titration.

Right flask: If further iodate solution were to be added after the endpoint (missed the endpoint).




15+


Part 2: Determination of Ascorbic acid Amount (%yield) in tablet solution *Procedure:*

1. Dissolve a single vitamin C tablet in 200 mL

of distilled water (in a volumetric flask if possible).

- 2. Take 20 mL of tablet solution and add 70 mL of distilled water in 250ml Erlenmeyer flask.
- 3. Add 10 mL of 1M H₂SO₄ and add 5 mL of starch (indicator).
- 4. Start titration with 0.05M lodine solution until you reach the end point.
- 5. Record the volume of iodine solution (titrant) used to reach the end point.
- 6. Calculate the (%assay) of ascorbic acid present in vitamin C tablet solution.





Part 1: Determination of Ascorbic Acid Amount (mg%) in Fruit Juice:

Notes:

- At the beginning, clean all glassware needed and label them clearly using permanent marker.
- Burettes are very delicate and expensive, make sure that you clean them while they are clamped to a stand.
- Rinse with distilled water, then rinse with the titrant solution.
- Do not pour any liquid inside the burette without using a liquid funnel and wearing safety goggles.



Ascorbic acid amount (mg //) in frint juice. Ex: the volume of (0.05 M) Iodine solution obtained was end }# of mol of AA = # of mol of I2] wint cup

point of AA / WWH AA = [MI2 * VI2] Exploises


weight of AA / WWH AA = [MI2 * VI2] Exploises ?? / 176.12 = 0.05 * 3 * 10⁻³ is latter weight of AA = 0.02649 ml = L End point 9-3 mg = paille 0.02649 * 0 = 26.4 mg in 10 ml of fruit Juice. mg /. -> X mg in 100 ml solution hint Comes

26.4 mg -> 10 ml procedure 2? -> 100 ml -> 264 mg/. in 200mles in collection in color of collection (in 200mles) in collection

طاف لماعانا البترية

The sequence of colors is clarified in the following figure:

- Left flask: before endpoint, iodine formed from added iodate reacts with ascorbic acid leaving the solution colorless.
- Centre flask: At the titration endpoint all the ascorbic acid has reacted, and the slight excess of added iodate reacts to give a pale blue color due to starch indicator. This is the endpoint of the titration.
- Right flask: If addition of iodate is continued after the endpoint, further starchiodine complex is formed, giving the solution a stronger blue-black color. (Missed the endpoint).



(Hint) Calculations:

The stoichiometric ratio of the reaction between ascorbic acid and iodine is 1:1.

Please refresh your memory with calculations of: # of moles, amount (mg%) and amount(%Assay) before attending the laboratory ession.

Part 3: Qualitative Determination of Ascorbic Acid:

molarity

Many chemical tests are performed to detect the presence of ascorbic acid.

الم مر موا

I. AqNO₃ test:

A.A "reducing agent" + AgNO3 HNO3 Ag ppt. (Grey ppt.) + AgNO3 + dehydroascorbic acid agent testular ste ació un, acid base Usli Not redox. oxidizing

II. NaHCO3 test:

A.A "acid" + NaHCO3 "base" -> CO2 + H2O + sodium ascorbate (Bubbling)

A.A + Benedict's reagent (Cu-2-citrate) heat Cu2O" brick red to brown ppt." + dehydroascorbic acid

الله Benedict's test:

A.A + Report بلدباته ماطع

Ascorbic Acid Amount (mg%) in Fruit Juice

Ex: the volume of (0.05M) Iodine solution obtained was 3 ml

Mol = weight / MWT

of mol of AA = # of mol of
$$I_2$$

Weight of AA / MWT AA = $M_{I_2} \times V_{I_2} \rightarrow E.P$
 $??/176.12 = 0.05 \times 3 \times 10^{-3}$
Weight of AA = 0.0264 g

 $0.0264 \text{ g} \times 10^{+3} = 26.4 \text{ mg}$ in 10 ml of fruit juice

Molecular weight for Ascorbic Acid= 176.12g/mol

mg% → X mg in 100ml solution

$$26.4 \text{mg} \rightarrow 10 \text{ ml}$$

?? $\rightarrow 100 \text{ ml} \rightarrow 264 \text{ mg}\%$

Ascorbic Acid Amount (%Assay) in tablet solution

Ex: the volume of (0.05M) Iodine solution obtained was 11.2 ml

Mol = weight / MWT

% Assay = Actual Amount Theoretical Amount

x 100%

One tablet contains 1 g of AA \rightarrow Theo. Amount 0.986 g \rightarrow Act. Amount

% Assay =
$$\frac{0.986}{19} \times 100\%$$