

EXTRACTION

A Separation and Isolation Technique

INTRODUCTION

مزيل مركب معين من مزيج
استخدم يفضل المادة يلي
Extraction تسمى
Solvent
بري افضل مادة صلبة
اصب مادة يلي بري استعملها
في الفصل تكون سائلة

Extraction is the separation of a substance from a mixture by means of a solvent that preferentially dissolves that substance. If the substance is extracted from a solid phase, the process is called solid-liquid extraction, as in the isolation of caffeine from tea leaves by means of hot water.

Extraction of a substance from a liquid phase is called liquid-liquid extraction. The most common applications of this latter technique are:

- The recovery of an organic product from a reaction mixture containing excess unreacted materials and by-products.
- Isolation of an organic substance from its natural source, such as a plant.

Liquid-Liquid Extraction. This is the most common type of extraction. It involves shaking the liquid mixture with an immiscible solvent which preferentially dissolves the desired compound. On standing, the two immiscible phases (usually organic and aqueous) form two separate layers (upper and lower) that can be separated by means of a separatory funnel.

hydrophilic
low organic
high polar

hydrophobic
High organic
low polar

The various solutes in the mixture distribute themselves between organic and aqueous phases according to their relative solubilities in each solvent. At equilibrium, the ratio of the concentration (C) or solubility of the substance in the organic phase, (C_o or S_o) to that in the aqueous phase (C_w or S_w) is called the distribution coefficient (K_D).

$$K_D = \frac{C_o}{C_w} = \frac{S_o}{S_w} \frac{\text{organic}}{\text{water}}$$

K_D كـ hydrophilic
Organic
via
low polar

A large distribution coefficient implies that the compound is much more soluble in the organic phase than in the aqueous one and, in this case, a single extraction suffices to remove the desired compound from the mixture. When K_D is small, it means that the compound distributes

K_D منخفض
↑ high polar hydrophilic
↓ low organic

كمية قريبة
للوامد 1.7
دليل تجريبي
هدول
سبة متساوية
في المركب، رات ترفع
بمركب الاكسام
المركب اكثر منه مرة
والمادة اكثر منه مرة

50
اذا كانت Polarity عالية او قليلة
راي يكون سهل جدا للمركب
من المركب لا يرفع
downlayer / upper layer

* إذا كاننا قريبي من الواحد لاني المفصل
 * انما جزء من المادة يعني اذمها الى اجزاء دوالي 30 دل 1 فصل اذم
 10ml 10ml 10ml ← كل ما انفصل على ما كانت عليه
 الفصل اذم .

more evenly in both phases, so that repeated extractions are required to recover such a compound from the aqueous mixture.

In general, it is more efficient to divide the total volume of extracting solvent over several extractions than to use the whole volume in a single extraction.

أقل من واحد / أكثر من واحد
 بكثير / بكثير

GENERAL EXPERIMENTAL CONSIDERATIONS

Choice of Solvent. A solvent used for extraction should have the following characteristics:

1. Immiscible with the liquid in which the solute is present.
2. Readily dissolve the solute to be extracted.
3. Extract little or none of the impurities and other compounds present in the mixture.
4. Non-flammable, non toxic, cheap and easily removable from the solute after extraction (*i.e.*, volatile).

المحلول متجانس
 benzene
 Dichloromethane
 غير ذائب في ماء
 غير سام
 غير قابل للاشتعال
 سهل التبخر
 غير قابل للاشتعال
 غير قابل للاشتعال
 سامة قليلة
 سهل التبخر

Salting-out. Extraction of organic compounds from aqueous mixtures is usually improved by saturating the aqueous phase with a salt such as NaCl or Na₂CO₃. This phenomenon is called salting-out and has the following effects:

1. Decreases the solubility of organic compounds in the saturated aqueous phase.
2. Decreases the solubility of the organic and aqueous phases in each other, thus improving their separation. This is particularly useful in breaking up emulsions.

Emulsions. In certain cases, the two immiscible phases do not separate cleanly into two distinct layers; instead, they form an emulsion which, once formed, is usually difficult to break. It is therefore advisable to prevent the formation of emulsions during extraction. This is best achieved by avoiding vigorous shaking of the layers whenever an emulsion is expected to form

* أغلب ما ينفصل ما بعض يتكد واضح
 51

* صعب التخلص .

* هي ما يكونه عندي استا و تزيد تليق 54 .

① عسيه يوع /

shake gently 51 قوي ليين جيد .

(e.g., when alkaline aqueous solutions are extracted with chloroform or dichloromethane). If an emulsion still forms one can often break it by:

1. Stirring the emulsified layer gently with a glass rod.
2. Saturating the aqueous layer with a salt.
3. Centrifugation.

Drying Agents. The organic phase often shows turbidity due to the presence of traces of water from the aqueous phase. Anhydrous CaCl_2 , MgSO_4 , or Na_2SO_4 may be used as drying agents which absorb the traces of water present in the organic phase. When dry, the organic phase becomes clear.

Acid-Base Extraction



Mixtures of organic acids and bases are commonly separated by acid-base extraction. Such compounds are converted to their salts by treatment with acid or base. Unlike the original compounds, the corresponding salts are usually soluble in water, thus enabling their transfer from the organic phase to the aqueous layer. After separation of the layers, the organic acid or base is recovered by neutralization of the aqueous layer. Since the acid or base is insoluble in water it precipitates out, and is collected by filtration. A flow diagram for such a separation is shown in Figure 16.

صلى انقش
منه emulsion
① ملع كبريتات الصوديوم
② مقرب في خفيته
ايضا بعض المواد

أخذ البقية
في لثة
التي يكون عندي
بعض قطرات
الماء معلقة
* الجيد drying agent
صلى ينقلها في

وزنه beaker
وهي بانه
راع بعض عندي
di chloromethane
traiffeme.
↑
Beaker

المزج وزنه
Beaker - Beaker




EXPERIMENTAL

MATERIALS NEEDED	
	Glassware: Beaker (600 mL), 2 beakers (100 mL), 2 Erlenmeyer flasks (100 mL), stand, wire gauze, clamp, ring, clamp holder, graduated cylinder (10 mL), separatory funnel (100 mL), cheesecloth 20x20 cm.
	Chemicals: 15 g Tea leaves, 10 g sodium carbonate, 30 mL dichloromethane, 1.0 g anhydrous sodium sulfate, 1 g benzoic acid, 1 g p-dichlorobenzene, 25 mL ether, 40 mL of 10% sodium hydroxide, 20 mL concentrated hydrochloric acid, anhydrous calcium chloride, blue litmus paper.

EXTRACTION OF CAFFEINE FROM TEA LEAVES

Caffeine is an organic compound present in the fruit and bark of some plants, as well as in tea leaves, coffee, cocoa and cola beans. The caffeine content in dried tea leaves is about 3-4%.



alkaloids
عباره عنه
نيتروجينيه
cyclic

Caffeine belongs to a family of basic, nitrogen-containing, cyclic compounds called alkaloids. It is a mild stimulant and is used as such in many drugs and analgesics. The solubility of caffeine at room temperature is 2.2 g/100 mL of water and 18 g/100 mL of chloroform.

In this experiment, you will extract caffeine from tea leaves with hot water. This treatment also extracts tannins, a class of acidic organic compounds, also present in the leaves. Sodium carbonate is used to remove the acidic tannins by converting them to water-soluble salts.

Procedure. (1) You will be provided with a large tea bag containing about 15 g dry tea leaves. (2) Place the tea bag in a 600 mL beaker, add 10 g of sodium carbonate and 150 mL of water, and boil the mixture



gently for 20 minutes. ^③ Cool the dark brown aqueous solution to room temperature and squeeze the tea bag to extract the liquid fully before discarding the bag.



- ④ Transfer the dark solution to a separatory funnel and extract twice with 15 mL portions of dichloromethane. Avoid vigorous shaking of ^⑤ the funnel since emulsions may form readily; instead, swirl the funnel gently or turn it upside down several times. [After each extraction ^⑥ drain the denser dichloromethane layer into a small flask.] [Dry the ^⑦ combined organic extracts with anhydrous sodium sulfate until the ^⑧ solution is clear. Decant the dichloromethane into a small beaker and evaporate to dryness over a water bath in the fume hood. Do not heat the residue any longer than necessary since caffeine decomposes readily.] [Weigh the crude caffeine and calculate its percentage in the tea leaves.] [Determine the melting point of your product.] ^⑩

① 150 mL water in 600 mL beaker

② Tea bag 15 g نبتة الأكلونديا  tea bag

③ 10 g sodium carbonate (10 غم) نسبة

④ Boil mixture gently for 20 minutes. انسود

⑤ cool the dark brown aqueous solution to room temperature.

⑥ Clamp holder + separatory funnel نصب

⑦ Transfer the dark solution into separatory funnel

→ Don't forget: to squeeze tea bag: to extract the organic layer fully before discarding it



⑧ Extract the solution twice 15 mL portions of dichloro methane

⑨ 15 mL dichloromethane. 55

* تصليح حقيفة

⑨ Drain the dichloro layer into small flask -methane

⑩ Repeat the extraction of solution with another 15 mL dichloro methane