EXPERIMENT 4

QUALITATIVE DETERMINATION OF PROTEINS AND AMINO ACIDS

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PROTEINS

- Proteins are one of the main constituents of the cell. Cells contain more than 3000 different proteins that play a variety of structural and functional roles in the cell.
- They are high molecular weight biopolymers of amino acids.
- The linkage that bind these amino acids is called a peptide bond and is formed by the condensation of the α -carboxyl group of one amino acid and an α -amino group of another

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DETECTION OF PROTEINS:


- Detection based on the following: Presence of peptide bond: Biuret reaction The nature of protein (coagulation by heating, precipitation by strong acid and heavy metals) Presence of certain amino acids
- Color reactions of amino acids
- Biuret test: Compounds containing one or two peptide bonds give a characteristic purple color when treated with dilute copper sulphate in alkaline solution. The positive reaction of Biuret reagent is formed by a complex of copper atom and four nitrogen atoms, two from each of two peptide chains.

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
Procedure:

- Place 15 drops of protein test solution in a test tube and 5 drops of 20% NaOH and mix. Add 2 drops of 0.2% CuSO₄ solution, warm if necessary Repeat with any amino acid instead of protein and write your notes.

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HEAT DENATURATION:

- Proteins are unequally stable when heated in aqueous solution, fractionation can be achieved by controlling heating.
- The presence of a substrate, a product, or an inhibitor often stabilizes a protein against heat denaturation.
- Denaturation consists of the unfolding polypeptide chain and loss of the compact structure. It is usually irreversible and can result from the application of heat, extremes of pH or the action of detergent.
- Denaturation therefore frequently causes precipitation of the protein
- Procedure:
 - Take 5 ml of protein solution in a test tube, heat to boiling, repeat using any amino acid instead of protein. What do you observe

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AMINO ACIDS:

- Amino acids are the main building blocks of proteins.
- They are relatively small molecules that characterized by presence of α -amino group, carboxylic group and a unique R group.
- Although they share the general structure, they differ in the polarity, size, solubility and electrical charge.

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NINHYDRIN TEST:

- Ninhydrin is a powerful oxidizing agent, it reacts with α -amino acid to decarboxylate the amino acid yielding an intensely colored blue-purple product, CO₂, water and the corresponding aldehydes.
- Ninhydrin yields a similar blue-purple reaction upon reaction with primary amines and ammonia. The reagent also reacts with amino acids such as proline to yield a yellow to red product. Under controlled conditions the color development is quantitative. Several other color reactions are generally specific to individual amino acids.
- These are particularly useful for locating specific amino acids, or peptides containing these amino acids on paper chromatograms. Amino acid residues for which good color tests have been achieved are arginine, cysteine, proline, hydroxyproline, tryptophan and tyrosine.

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
Arginine, Cysteine, Proline, Hydroxyproline, tryptophan, tyrosine

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NINHYDRIN TEST:


- **Procedure:**

- In different test tubes containing 15 drops different amino acids and protein, add 0.5 ml of 0.1% ninhydrin, heat to boil for 1-2 minutes, allow cooling Repeat with ammonia instead of protein.

MILLON'S REACTION:

- Specific for tyrosine and proteins containing it.
- Millon's reagent is a mixture of mercurous and mercuric nitrates in nitric acid which reacts with the phenolic ring of tyrosine to form red mercury phenolate.
- Inorganic salts in large amounts interfere with the reagent by precipitating mercury.

- Procedure:
- In a test tube, add 6-8 drops of millon's reagent to 2 ml of tyrosine solution. A white precipitate may form ($\text{HNO}_3\text{Hg}^{++}$) Heat the tube to the boiling point, if a brick red color developed add 2-3 drops of the reagent and heat again. Write your observation

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
HOPKINS-COLE REACTION:

- The test is specific for proteins containing tryptophan due to the interaction of an aldehyde of glyoxylic acid with indole nucleus.
- Glyoxylic acid is prepared from glacial acetic acid by being exposed to sunlight

- Procedure:
- Mix 3 ml of tryptophan solution with an equal volume of Hopkins-Cole reagent in a test tube, gently add 3-5 ml of conc H_2SO_4 along the side of the tube (DON'T MIX) Notes the result

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
SULFUR REACTION:

- Sulfur containing amino acids on heating with KOH, a black lead sulphide precipitate is formed by reacting with lead acetate

(Detection of sulphur containing amino acids)

- Procedure:**
- In two separate test tubes, place 2ml of cysteine solution, add 2ml of 10% KOH and mix. Add 2 drops of 1% lead acetate solution, boil the mixture of the tube. The color will deepen into black if sufficient sulfur is present.

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
NITROPRUSSIDE TEST:

- It is specific for cysteine, the only amino acid contains sulfhydryl group (SH) which reacts with sodium nitroprusside in the presence of ammonia to give a purple color

Nitroprusside Test (Detection of Cysteine)

- Procedure:
- In two separate test tubes, place 2 ml of solution (cysteine or cystine). Add 2ml of 20% NaOH Add 5-6 drops of 2% Na nitroprusside solution What do you observe

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