# EXPERIMENT 10 GLUCOSE 1- DEHYDROGENASE ACTIVITY ASSAY

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
\* (( G6P )) dehydrogenase ... damage Red Blood Cell من الإبريق

\* Hexokinase: an enzyme that phosphorylates.  
\* ATP      \* NAD<sup>+</sup>

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1 vial  
\* GDH positive control (lyophilized)  
\* NADH standard  $\rightarrow$  (lyophilized)

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\* ② pipetted up & down Several time.  $\rightarrow$  \* حتى powder يذوب \*  
.. Mixing & Shaking \*  
مربى

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\* ② Aliquot & Store at -20°C.

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\* the GDH positive Control  $\rightarrow$  less stability  $\rightarrow$  2 months at -20°C after reconstitution ..

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Reconstituted  $\rightarrow$  NADH (10 mM)  
 $\rightarrow$  Supplied Glucose (2M)  $\rightarrow$  Stable for up 6 months at -20°C.

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... Saia dwaitkat

# GENERAL TROUBLESHOOTING GUIDE:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> <li>• Use of ice cold assay buffer</li> <li>• Omission of a step in the protocol</li> <li>• Plate read at incorrect wavelength</li> <li>• Use of a different 96-well plate</li> </ul>	<ul style="list-style-type: none"> <li>• Assay buffer must be at room temperature</li> <li>• Refer and follow the data sheet precisely</li> <li>• Check the wavelength in the data sheet and the filter settings of the instrument</li> <li>• Fluorescence: Black plates (clear bottom), Luminescence: White plates, colorimetric: clear plates</li> </ul>
Samples with erratic readings	<ul style="list-style-type: none"> <li>• Use of incompatible sample type</li> <li>• Samples prepared in different buffer</li> <li>• Samples were not deproteinized (if indicated in data sheet)</li> <li>• Cell/tissue sample were not completely homogenized</li> <li>• Samples used after multiple freeze-thaw cycles</li> <li>• Presence of interfering substance in the sample</li> <li>• Use of old of inappropriately stored samples</li> </ul>	<ul style="list-style-type: none"> <li>• Refer data sheet for details about incompatible samples</li> <li>• Use the assay buffer provided in the kit or refer data sheet for instruction</li> <li>• Use the 10-kDa spin cut-off filter or PCA precipitation as indicated</li> <li>• Use the Dounce homogenizer (increase the number of strokes), observe for lysis under microscope</li> <li>• Aliquot and freeze samples if needed to use multiple times</li> <li>• Troubleshoot if needed</li> <li>• Use fresh samples or store at correct temperatures until use</li> </ul>

# GENERAL TROUBLESHOOTING

Lower/higher readings in samples and standards	<ul style="list-style-type: none"><li>• Improperly thawed components</li><li>• Use of expired kit or improperly stored reagents</li><li>• Allowing the reagent to sit for extended times on ice</li><li>• Incorrect incubation times or temperature</li><li>• Incorrect volumes used</li></ul>	<ul style="list-style-type: none"><li>• Thaw all components completely and mix gently before use</li><li>• Always check the expiry date and store the component appropriately</li><li>• Always thaw and prepare fresh reaction mix before use</li><li>• Refer data sheet and verify correct incubation times and temperature</li><li>• Use calibrated pipettes and aliquot correctly</li></ul>
--	--	--

# GENERAL TROUBLESHOOTING GUIDE:

<p>Readings do not follow a linear pattern for standard curve</p>	<ul style="list-style-type: none"> <li>• Use of partially thawed components</li> <li>• Pipetting error in the standard</li> <li>• Pipetting error in the reaction mix</li> <li>• Standard stock is at an incorrect concentration</li> <li>• Calculation error</li> <li>• Substitution reagents from older kits/lots</li> </ul>	<ul style="list-style-type: none"> <li>• Thaw and resuspend all components before preparing the reaction mix</li> <li>• Avoid pipetting small volumes</li> <li>• Prepare a master reaction mix whenever possible</li> <li>• Pipette gently against the wall of the tubes</li> <li>• Always refer the dilutions in the data sheet</li> <li>• Recheck calculations after referring data sheet</li> <li>• Use fresh components from the same kit</li> </ul>
<p>Unanticipated results</p>	<ul style="list-style-type: none"> <li>• Measured at incorrect wavelength</li> <li>• Samples contain interfering substances</li> <li>• Use of incompatible sample type</li> <li>• Sample readings above/below the linear range</li> </ul>	<ul style="list-style-type: none"> <li>• Check the equipment and filter setting</li> <li>• Troubleshoot if it interferes with the kit</li> <li>• Refer data sheet to check if the sample is compatible with the kit or optimization is needed</li> <li>• Concentrate/dilute sample so as to be in the linear range</li> </ul>