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
Experiment 4 Qualitative Determination of Protein and amino acids

Protein is an important macronutrient essential for survival. They are an important constituent of cells and hence are present in all living bodies. Cells contain more than 3000 different proteins that play a variety of structural and functional roles in the cell.

10-35% of calories should come from protein. Protein is found in meats, poultry, fish, meat substitutes, cheeses, milk etc. They contain carbon, hydrogen, oxygen, nitrogen, phosphorus and Sulphur in some cases. Proteins are large biological molecules of alpha-amino acids.

Amino acids are the main building blocks of proteins. They are relatively small molecules that characterized by presence of alpha-amino group, carboxylic group and a unique R group. the amino group is attached to alpha carbon, which are crystalline in nature and exist as zwitterions (it has both positive and negative)

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Although amino acids share the general structure, they differ in the polarity, size solubility and electrical charge.

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
Amino acids form a specific type of linkage known as **peptide linkage**, amino acid molecules undergo a condensation reaction. the condensation of the alpha-carboxyl group of one amino acid and an alpha-amino group of another. The products formed are classified as; depending on the number of amino acid molecules involved in the condensation reaction.

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Proteins differ primarily from each other in their amino acid sequence. There are about 20+ amino acids here. Some amino acids are not produced by the body and are delivered by diet. They are called amino acids, which are **essential**.

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
This test is used to detect the presence of peptide bond (not less than two peptides) .

When treated with copper sulphate solution in presence of alkali (NaOH or KOH), protein reacts with copper (II) ions to form a violet coloured complex called biuret.

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B. Heat Denaturation:

Proteins are unequally stable when heated in aqueous solution, fractionation can be achieved by controlling heating. The presence of a substrate, a product, or an inhibitor often stabilizes a protein against heat denaturation. Denaturation consists of the unfolding polypeptide chain and loss of the compact structure. It is usually irreversible and can result from the application of heat, extremes of pH or the action of detergent. Denaturation therefore frequently causes precipitation of the protein.

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
Ninhydrin is a powerful oxidizing agent, it reacts with **α -amino acid** to decarboxylate the amino acid yielding an intensely colored blue-purple product, **CO₂**, **water** and the **corresponding aldehydes**.

Ninhydrin Test

© Bivus.com

Ninhydrin yields a **similar blue-purple reaction** upon reaction with **primary amines and ammonia**. The reagent also reacts with amino acids such as **proline to yield a yellow to red product**.

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This test give a **positive result** by **only amino acids and proteins** which **contain free – NH₂ groups (α amino group)** in their structure. Amino acid residues for which good color tests have been achieved are **arginine, cysteine, proline, hydroxyproline, tryptophan and tyrosine**. The colour intensity is **proportional to the concentration of amino acid i.e. ammonia from the amino groups**.

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
B. Millon's Test:

Millon's reagent is a **concentrated HNO₃** reagent that **dissolves mercury**.

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The test by Millon is specific to **phenol** that contains structures. Specific for **tyrosine** (tyrosine is the only common phenolic amino acid) and **proteins containing it**.

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A red precipitate or a red solution is regarded as a **positive test** as a consequence of the reaction. **Inorganic salts in large amounts interfere** with the reagent by **precipitating mercury**. A **yellow HgO precipitate is NOT a positive response**, but it usually shows that the solution is **too alkaline**.

This test can be applied to **tyrosine, phenylalanine**.

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
C. **Hopkin's Cole Test:**

Hopkin's Cole test is a specific test used for the detection of **indole ring** and thus, **tryptophan** in proteins. The test is also termed as '**glyoxylic acid test**' as the reagent contains **glyoxylic**. Glyoxylic acid is prepared from **glacial acetic acid** by being **exposed to sunlight**.

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In the presence of **concentrated H2SO4**, the **indole tryptophan group** reacts with **glyoxylic acid** to give a **purple color**. The **H2SO4** added to the reagent helps to **stabilize the glyoxylic acid** and **prevent its decomposition** and the release of **carbon dioxide**.

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
Sulfur containing amino acids containing sulfur, S-S group in cysteine, and S-H group like cysteine and cystine on heating with **KOH**, a **black lead sulphide precipitate** is formed by reacting with lead **acetate**

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E. Nitroprusside's test

The nitroprusside's test is specific for cysteine, the only amino acid containing a sulfhydryl group (—SH). The group reacts with nitroprusside in an alkaline solution to yield a red complex. a free sulfhydryl group (-SH) that is able to react with the nitroprusside in the presence of excess ammonia (NH4OH).

Reactions:


Fig: Na-nitroprusside reaction with Cys amino acid.

Practical work

A. Biuret test:

1. Place 15 drops of each sample (unknown) test solution in a separate test tube.
2. add 5 drops of 20% NaOH and mix.
3. Add 2 ml of 0.2% CuSO₄ solution, **warm if necessary**

Biuret Test

Negative biuret test

positive biuret test (polypeptide chain chelates with a copper ion)

B. Heat denaturation:


Take 5 ml of each sample (unknown) test solution in a separate test tube, heat to **boiling**.

What do you observe?

C. Ninhydrin Test:

1. Take 2 ml of each sample solution in a separate test tube.
2. Add 0.5 ml of 0.1 % Ninhydrin reagent the test tubes.
3. Place the test tubes in the water bath for 5 minutes and then allow cooling to room temperature.
4. Observe the formation of a deep blue/purple colour indicates the presence of amino acids.


Result:

- ✓ blue-purple reaction
- ✓ yellow or yellow to red color

Interpretation

1. This reaction suggests that amino acids, other amines, and ammonia are present in the test material.
2. Negative: not a.a , not protein
3. Positive purple: aa or protein with free amino gp
4. Positive yellow: proline aa or hydroxyproline.

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
D. Millon's test

1. Take 2 ml of each sample solution or unknown solution is taken in a separate test tube.
2. To this, about 2 ml of Millon's reagent is added. Observe the formation of red brick precipitate
3. if red colored precipitate is not observed immediately. the test tubes are then kept in the water bath for about 2 minutes The tubes are then observed for the formation of the colored precipitate.

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
Positive result: A positive result in the Millon's test is demonstrated by the formation of a red brick or pink colored precipitate. This indicates the presence of tyrosine or tyrosine containing protein.

Negative result: A negative result in the Millon's test is demonstrated by the absence of colored precipitate in the test tube. This indicates the absence of tyrosine or tyrosine-containing protein.

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E. Hopkins-Cole test

1. In two separate test tubes, take 2 ml of Hopkins-Cole reagent and 2 ml of each sample liquid are taken.
2. To this, concentrated H_2SO_4 is added along the sides of the test tube held at a slanting position. Two distinct layers of liquid are to be formed without mixing. (Note: Mouth of the tube must point away from the face, as we should be careful while adding the sulfuric acid.)
3. The test tube should be observed for the formation of a purple colored ring at the interface of two layers.

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Interpretation of result

Positive result: A positive result is represented by the formation of a purple-colored ring at the junction of two layers. This indicates the presence of tryptophan-containing proteins.

Negative result: A negative result is represented by the absence of a purple-colored ring in the test tube. This indicates the absence of tryptophan-containing proteins.

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F. Lead sulfur test

1. In two separate test tubes, place 2 ml of sample solution add 2 ml of 20% NaOH and boil for a minute.
2. Cool it and add a drop of 20 % lead acetate solution.
3. Observe the formation of black lead sulfide precipitate.

Note: Carry out this test in the hood if possible.

Positive test: Formation of black precipitate indicate the presence of sulfur-containing amino acid.

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
Negative test: No formation of black precipitate indicates the absence of sulfur-containing amino acid.

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G. Nitroprusside test:

1. In two separate test tubes, place 2 ml of sample solution.
2. Add 0.5 mL nitroprusside solution and shake thoroughly.
3. Add 0.5 mL of 20% Sodium hydroxid .
4. Observe the color change.

Result interpretation

Positive test: Formation of red color indicate presence of cysteine.

Negative test: No formation of red color indicates absence of cysteine.

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