




إعداد الصيدلاني/ة: محمد يأسين








Some active constituents need to be processed (slow drying is allowed)

Liberation of active vanillin by enzymatic process from Gentiana Luten and



3. Acelleular products

They are plant isolates like volatile oil, balsams, resins and fixed oils.

Methods to extract acellular products (Volatile oils):

- 1. destructive distillation: wood of Scotch pine (used for dermititis)
- 2. Steam distillation using Clavenger tube: plant material is heated with water and steam is passed through the mixture and condensed. This needs a round-bottomed flask لله that contains water or alcohol (usually water), the plant is inserted inside it with the water, then it's attached into a condenser, also a collecting tube should be there to collect the extracted oil, distillation also requires a heat source, if the treated material is very sensitive for heat (thermo labile) the steam distillation can't be used then, but other methods like:
- ~ Hydrosteam distillation: here the plant is separated from the water by a membrane, the plant part is inserted to the membrane and by water heating, the vapor passes this membrane to the plant material, the v.o. is then extracted.

3. Acelleular products

- 3. Expression العصر used pharmaceutically for citrus oils only, fruits are turned by a stainless steel perforated roll, each pore is attached into an opening ending into a collector, upon turning these pores are induced.
- 4. Solvent extraction, using alcohol, not used in pharmaceutical industry.
- 5. Enzymatic hydrolysis: very uncommon, some V.O. that has a sugar part force to use this method to extract the non-sugar part (previously mentioned example: Vanilla).
 - (previously mentioned example: Vanilla).

 6. Enfleurage: in expensive perfume industry, an old method that uses sheets of waxes or fatty solid materials, the plant is inserted into these glassy sheets, then put in woody frames for a week or 2 weeks in a dark area, then the plant part used is removed and the sheets then are treated with alcohol to produce a concentrate of these v.o.








tside cell

3. Acelleular products VOLATILE OIL STEAM DISTILLATION



Clavenger tube

Figure 1. Clevenger-type apparatus — Appareil de type Clevenger.

3. Acelleular products VOLATILE OIL STEAM DISTILLATION

Hydro-steam distillation

Clavenger tube




Figure 1. Clevenger-type apparatus — Appareil de type Clevenger.

عملية التقطير لتحضير الزيوت الطيّارة VOLATILE OIL STEAM DISTILLATION

- 2/2/20
- Principle: When water is mixed with another immiscible liquid and both are boiled together boiling will occur at lower temperature than the water's boiling point
- Plant material is heated with water and steam is passed through the mixture and condensed were is adifferent
- Disadvantage of steam distillation:
- 1. unstable material may be destroyed
- 2. other substances in the plant may affect the composition of the oil
- Ex. In Citrus spp. The acidic nature cause a number of acidic nature cause and acidic nature cause acidic nature cause and acidic nature cause acidic nature cause acidic natu

white reposition of the rest of the reposition of the reposition of the collect of the rest of the collect of the rest of the

4. Galenical preparations and extracts

- Galenical: مستحضر طبي من اصل نباتي
- Galenical preparations: Are extracts of plant materials used in therapy.
- Methods of extraction to any plant:
- A. Infusion: drug stood in hot or cold water or other solvent for a short period of time.
 - Ex. Preparation of plants containing volatile oil.

(decation), shall being we led the less him is the little of the little

الما بالنقع (infusion) الزيو تح الطيارة تبعر بالماء.


Methods of extraction to any plant:

colle de

• **B. Decoction**: drug boiled with water or solvent.

e Ex.: roots, rhizomes and barks

C. Maceration: prolonged infusion (2-3 weeks), often in aqueous alcohol



• E. Hot extraction: repeated infusion with hot solvent example: Soxhlet extraction



like offer machine

Citter water


Not water



> COMMON METHODS OF PLANT EXTRACTION

Percolation:

In this method, the plant material is subjected to a slow flow of fresh solvent at intervals until sufficient active ingredient is extracted.






E. Soxhlet extraction:

- Soxhlet extraction is only required where the desired compound has a limited solubility in a solvent, and the impurity is insoluble in that solvent.
 - If the desired compound has a high solubility in a solvent then a simple filtration can be used to separate the compound from the insoluble substance.
 - The advantage of this system is that instead of many portions of warm solvent being passed through the sample, just one batch of solvent is recycled.
 - This method cannot be used for thermolabile compounds as prolonged heating may lead to degradation of compounds

Extraction: Withd of soluble constituents from

Extraction: Withdrawal treatment with solvent



Concentration of filtrate to dryness


عارة لزج من الحول

- The filtered plant part (liquid portion) must be reduced to paste-like, cakey, and finally powdered form through standard techniques.
- The essence of concentrating filtrate to dryness is to:
- Calculate the active ingredient present after the extraction by calculation of percentage yield.
- To obtain dose of plant extract in milligram or gram that can be further constituted in a known volume of vehicle for experimental animals dosing.

The enriched extract from maceration or extractors, concentrated by;

- Peding into a water bath, oven, hot plate, heating mantle, ...etc to produce a thick concentrated extract free from solvent.
- ② Or the concentration could be done using rotary evaporator in order to remove excess solvents from samples by applying heat to a rotating vessel at a reduced pressure.
- A freeze dryer or lyophilizer can also be used to concentrate and dry the filtrate extract (only aqueous) at reduced temperature

COMMON METHODS OF CONCENTRATING FILTRATE TO DRYNESS





Heating mantle



Rotatory evaporator

GRAPHICAL OUTLINE

